

EVERYTHING HAPPENS FOR A REASON...

So the saying goes and every year within the US 100,000 people die from DVT and yet Intermittent Pneumatic Compression (IPC) can reduce this risk by up to 60%... Then why are not more patients and clinicians using IPC.



.... and their reasons include

- Static connection to the bedside pump.
- Noise from every cycle reduces rest.
- Nurses are frustrated with excessive time spent, removing/reapplying limb sleeves.
- Hoses and limb sleeves slipping down may lead to accidental falls and ineffectiveness.
- Limited battery life only lasts for a few hours before recharging.
- High cost both financially and to the environment with the disposal of sleeves.
- Unknown patient compliance and device location resists use at home.



What is the reason behind Legacy

Only when new technology becomes available can change happen. Legacy uses a novel piezo pump system which provides compression to the calf mimicking ambulation, just like other Intermittent Pneumatic Compression devices. This technology however allows for the smallest, lightest, portable IPC, which operates in total silence. It also allows for extended battery life and remote connection allowing you to track the device and usage. We've also thought about the environment and limited the amount needed to be thrown away.

...And that is the reason why we developed Legacy

The Legacy IPC can be used straight from the box and includes:-


- 2 x Pumps
- 2 x Bladders (LH and RH)
- 2 Comfort Sleeves (LH and RH)
- 1x Power supply and Lead
- User Manual

All parts are reusable except the Comfort Sleeves which are single use



The Legacy can be connected via Bluetooth to any iOS or Android enabled device. Download the mobile App from the website: www.Legacy-IPC.com to monitor the status of the device.



Silent	1-2dBA over ambient	Only totally silent IPC on the market
Compact	4.17 x 1.77 x 0.67 inch	Smallest IPC on the market
Lightweight	0.154 lbs	Lightest IPC on the market
Long run time	>24 hours on a single charge	Only portable IPC which can run for more than 24 hours on a single charge.
Connectivity		Only IPC with remote connection to any iOS and Android device showing device status or later retrieval of information.
Effectiveness	45mmHg Compression	45mmHg ±5mmHg compression every minute
Reduced Costs	50% reduction	Running costs compared to the Flowtron and Kendall SCD
Environment	66% reduction	Reduction in plastic waste compared to Flowtron and Kendall SCD sleeves
Compliance	Home or Hospital	Can be used in both the hospital or community care



For ordering or for more information please contact:-



www.VasoCARE.com
 800-256-9979 (Corporate)
 800-455-5150 (Fax)
john@vasocare.com